aquapick®

"NAPREDNI ORALNI IRIGATOR"

Neophodan za oralno zdravlje
Novi Aquapick

aquapick

Made in Korea

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A woman with blonde hair is running outdoors. She is wearing a black sports bra and black leggings. The image is annotated with red dashed lines and labels in white rounded rectangles. The labels point to the brain, respiratory system, heart and blood vessels, and stomach. The background is a blurred green forest.

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- *Journal of American gerontology* (Denmark Research Team) -

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- See the February, 2002, issue of the *British Journal on the Periodontal Membrane* -

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- See the March issue of the *American Heart Association Journal on 'Circulation'* -

MATERICA

7,8 PUTA VEĆI RIZIK ZA PREVRREMEN POROĐAJ U ODNOSU
NA ZDRAVE LIUDE

- See data collected by the *International Dental Research Association* -

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- See data collected by the *Taiwanese Dental Health Association* -

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EU certification

German certification for meeting the European Safety Act

American safety certification

Korean certification for safe electric appliances

Japanese certification

European certification for meeting the approved electromagnetic standard

Recommended by the Swiss Agency for Health

Green Incentive Priority office

FDA certification

AQ-300
Made in Korea

MIKODENTAL DENTAL DEPOT

Glavno sedišta firma:
Stojana Novakovića 19
15000 Šabac/Srbija
Tel/fax: +381 15 345 123
Tel: +381 15 300 261

MIKODENTAL BG

Gandijeva 130a
11070 Novi Beograd/Srbija
Tel: +381 11 411 89 48
Tel: +381 64 6438154

MIKODENTAL BH

Derviša Numića 11
71 000 Sarajevo/BiH
Tel/fax: +387 33 712 111
Tel: +387 61 244 814

MIKODENTAL CROATIA

CIRKON-DENTAL
Smičklasova 8
47000 Karlovac
Tel: +385 91 5464696

aquapick
www.aquapick.com

Aquapick Co., Ltd.

C-1006, Woorim Lion's Valley, 168, Gasan digital 1-ro, Geumcheon-gu, Seoul, South Korea
Tel: +82,1588,2809 Fax: +82,2,6116,2800 E-mail: contact@aquapick.com